

STL Export Instructions for iTero Scanners

1. Log into www.myaligntech.com.
2. On the top left there is a box under **Quick Search**. Select the drop-down menu and choose **Patient Name**.
3. Type the patient name in and they will come up under **Orders**. Click on their **ID #** in blue and this will take you to the order information page.
4. **Click on Export (OrthoCAD 3.5 or higher)**. If prompted with a security warning, click **Allow**.

Click "Export", which requires OrthoCAD 3.5 or higher, to begin the export process.

Order Information

Order ID: 7007403 First Name: A * Company: QACompany - Orthodontic
Order Code: 7K6T8F Last Name: Demo * Doctor: Dr. QA Exclusively, QAPerso
Case Type: iCast Chart Number: Mode01 Doctor License: 111
Registration Date: 2011/08/14 Acquisition Date: 2011/08/14 Tracking Number:
State: Completed Order Status: Active
Notes:

Actions: [Open File](#) [Export \(OrthoCAD 3.5 or higher\)](#) [Change Order](#) [Image Management](#)

Order Details ID	Item	Quantity	Ship To	Bill To	Due Date	Shipping Date	Delivery Date	Status
7007404	iCast	1.00	QACompany - Orthodontic Office - iOC		2011/08/18	2011/08/17		Completed

5. Now a study model will pop up and a window that says **OrthoCAD Export**.
6. **Export type:** Choose the **Open Shell** option.


OrthoCAD Export

Export Type: Open Shell
Data Format: Hollow ABO Base
Folder Name: 7K6T8F


It is always up to the doctor to decide the treatment plan for the patient.
Third party providers may provide guidance on the export type that would be most compatible for your intended use.

Export

7. **Data Format:** Choose **Two Files (arches oriented in occlusion)**.




8. **Folder Name:** Change the number to the last name of your patient.



9. Click **Export**

10. A new window will open that says “Model data was exported to c:\cadent\export/(patients last name).” Click **Yes** to open the export folder.



11. A new window will open with several files in it. The first two files are in the .stl format that we use. The first one has a "l" for lower and the second "u" for upper and both are designated with a golden certificate. You can now upload those files to DPS. To upload to DPS, call us at 800-332-3341.

