






A Guide to Selecting Metal Free Restorations



CERTIFIED
DENTAL
LABORATORY



	Nautilus Full Contour 	Nautilus Core 	Nautilus Elite 	Simple Z 	e.max 
Material	Monolithic zirconia	Porcelain fused to zirconia	Monolithic zirconia	Monolithic zirconia	Lithium disilicate
Average Strength	1400 MPa	1400 MPa <i>(coping only)</i>	700 MPa	1400 MPa	400 MPa
Indications	Bicuspids Posterior Singles Posterior Bridges	Anterior Singles Bicuspids Bridges Posterior Singles	Anterior Singles Anterior Bridges Bicuspids	Posterior Singles <i>(primary anatomy only)</i>	Anterior Singles Anterior Bridges Bicuspids
Maximum Bridge Span	2 pontics <i>(recommended)</i>	2 pontics <i>(recommended)</i>	1 pontic	N/A	1 pontic
Shade Options	Vita Classic Vita 3D	Vita Classic Vita 3D Custom shading and characterization	Vita Classic Vita 3D Limited characterization	Vita Classic	Vita Classic Vita 3D Custom shading and characterization
Appropriate for Titanium Abutments	Yes	Yes	No	Yes	No
Cementation	Conventional Cementation	Conventional Cementation	Conventional Cementation	Conventional Cementation	Conventional Cementation or Bond

ADA Codes for Metal Free Restorations:

Single Crown D-2740

Bridge Abutment D-6740

Bridge Pontic D-6245

CERTIFIED
DENTAL
LABORATORY



DENTAL PROSTHETIC SERVICES, INC
1900 51st Street NE • Cedar Rapids, Iowa 52402
319-393-1990 • 800-332-3341 • www.DPSdental.com