

# *A Guide to Selecting Dentures*



## Step 1 - Select a Denture Line

	<b>Value Line</b> <i>Best for immediate placement* or for budget conscious patients.</i>	<b>Integrity Line</b> <i>A balance of function, fit, esthetics, and price.</i>	<b>Ultra Line</b> <i>Our most customized denture, offers the highest level of esthetics, fit, and function.</i>
<b>Material and Finish</b>	Ivocap® acrylic Smooth finish Smooth palate	Ivocap® acrylic Naturally contoured gingiva Anatomical palate	Ivocap® acrylic Naturally contoured gingiva Customized gingiva coloring available Anatomical palate
<b>Articulation</b>	Stratos 100	Stratos 100	Stratos 100
<b>Set Up</b>	Basic set up	Esthetic tooth arrangement	Custom anterior esthetic tooth arrangement
<b>Resets Included</b>	No	1 Reset - if necessary	2 Resets - if necessary
<b>Base and Bite Registration</b>	Light cure base, bite block	Your choice of base and bite registration	Your choice of base and bite registration
<b>Personal ID</b>	No	No	Yes

## Step 2 - Select a Tooth Type

	<b>Ivostar®/Gnathostar®</b> <i>For value-driven restorations.</i>	<b>BlueLine® DCL</b> <i>For any removable restorations.</i>	<b>Phonares® II</b> <i>For implant restorations or highly esthetic and high wear considerations.</i>
<b>Material</b>	PMMA	Double Cross Linked (DCL)	Nano Hybrid Composite (NHC)
<b>Shade</b>	16 Vita Classic Shades, Bleach BL2 and BL4	16 Vita Classic Shades, Bleach BL1 to BL4	16 Vita Classic Shades, Bleach BL1 to BL4
<b>Anterior Moulds</b>	15	42	18

### ADA Codes:

Complete Denture - Maxillary D5110  
 Immediate Denture - Maxillary D5130  
 Adjust Complete Denture - Maxillary D5410  
 Adjust Complete Denture - Mandibular D5411

Complete Denture - Mandibular D5120  
 Immediate Denture - Mandibular D5140

\*Immediate set up fee - \$3.00 per tooth



DENTAL PROSTHETIC SERVICES, INC  
 1900 51st Street NE  
 Cedar Rapids, Iowa 52402  
 319-393-1990 • 800-332-3341  
[www.DPSdental.com](http://www.DPSdental.com)

