

# *A Guide to Selecting Partial Dentures*



Step 1 - Select a Partial Type	Acrylic Partial	Flexible Partial	Cast Frame Partial	Ultaire™ AKP
<b>Material</b>	Wrought wire Ivocap® acrylic	Thermoplastic - Duraflex Nylon resin - Valplast	Chromium Cobalt Ivocap® acrylic	Polymer Ivocap® acrylic
<b>Articulation</b>	Hand or nonadjustable	Stratos 100 - adjustable	Stratos 100 - adjustable	Stratos 100 - adjustable
<b>Retention</b>	Tooth born	Tissue born	Tooth born	Tooth born
<b>Indications</b>	Temporary or provisional restorations Case plans involving future removal of teeth	Most partial applications Metal sensitivity	Most partial applications	Most partial applications Metal sensitivity
<b>Contraindications</b>	Long term use	Immediate applications Flat ridge Minimal vertical clearance Poor soft tissue support	Metal sensitivity Lack of suitable teeth in the arch Severe periodontal conditions	Lack of suitable teeth in the arch Severe periodontal conditions
<b>Clasping Color</b>	None or wrought wire	Gingival pink	Chromium Cobalt	Bone
<b>Add Teeth</b>	Yes	Yes	Yes	Not at this time

Step 2 - Select a Tooth Type	Ivostar®/Gnathostar® <i>For value-driven restorations.</i>	BlueLine® DCL <i>For any removable restorations.</i>	Phonares® II <i>For implant restorations or highly esthetic and high wear considerations.</i>
<b>Material</b>	PMMA	Double Cross Linked (DCL)	Nano Hybrid Composite (NHC)
<b>Shade</b>	16 Vita Classic Shades, Bleach BL2 and BL4	16 Vita Classic Shades, Bleach BL1 to BL4	16 Vita Classic Shades, Bleach BL1 to BL4
<b>Anterior Moulds</b>	15	42	18

### ADA Codes:

Cast Metal Partial Denture - Maxillary D5213

Resin Base Partial Denture - Maxillary D5211

Flexible Partial Denture - Maxillary D5225

Unspecified Removable Prosthodontic Procedure, By Report - D5899 (Ultaire)

Cast Metal Partial Denture - Mandibular D5214

Resin Base Partial Denture - Mandibular D5212

Flexible Partial Denture - Mandibular D5226



DENTAL PROSTHETIC SERVICES, INC  
1900 51st Street NE  
Cedar Rapids, Iowa 52402  
319-393-1990 • 800-332-3341  
www.DPSdental.com

