

# Abutment Torque Reference Sheet



Manufacturer	Implant System	Torque (Ncm)		Driver Tip
<b>BioHorizons</b>				
	Internal/Tapered Internal 3.5, 4.5, 5.7	30 Ncm	⬡	0.050" (1.25 mm) Hex
	External 3.5, 4.0, 5.0	30 Ncm	⬡	0.050" (1.25 mm) Hex
<b>Dentsply Sirona</b>				
	Ankylos® C/X 3.5, 4.5, 5.5, 7.0 (indexed or non-indexed)	15 Ncm	⬡	1.0 mm Hex
	Astra Osseospeed™ EV 3.0 Green, 3.6 Purple, 4.2 Yellow, 4.8 Blue, 5.4 Brown	25 Ncm	⬡	0.050" (1.25 mm) Hex
	Astra Osseospeed™ Profile EV 4.2 Yellow, 4.8 Blue	25 Ncm	⬡	0.050" (1.25 mm) Hex
	Astra Osseospeed™/Osseospeed™ TX 3.0 Yellow	15 Ncm	⬡	0.050" (1.25 mm) Hex
	Astra Osseospeed™/Osseospeed™ TX 3.5, 4.0 Aqua	20 Ncm	⬡	0.050" (1.25 mm) Hex
	Astra Osseospeed™/Osseospeed™ TX 4.5, 5.0 Lilac	25 Ncm	⬡	0.050" (1.25 mm) Hex
	Astra Osseospeed™ 4.0 Blue (Previous)	25 Ncm	⬡	0.050" (1.25 mm) Hex
	Astra Osseospeed™ TX Profile 4.5, 5.0	25 Ncm	⬡	0.050" (1.25 mm) Hex
	XIVE® S 3.0, 3.4, 3.8, 4.5, 5.5	24 Ncm	⬡	1.22 mm Hex
	FRIALIT™ 3.4, 3.8, 4.5, 5.5	24 Ncm	⬡	1.22 mm Hex
<b>Nobel Biocare®</b>				
	Angulated Screw Channel	35 Ncm	⬡	Nobel Biocare® Omnigrip
	NobelActive™ 3.0	15 Ncm	☆	Nobel Biocare® Unigrip
	NobelActive™/NobelReplace™ Conical Connection NP 3.5	35 Ncm	☆	Nobel Biocare® Unigrip
	NobelActive™/NobelReplace™ Conical Connection RP 4.3, 5.0	35 Ncm	☆	Nobel Biocare® Unigrip
	NobelReplace™ NP 3.5	35 Ncm	☆	Nobel Biocare® Unigrip
	NobelReplace™ RP 4.3	35 Ncm	☆	Nobel Biocare® Unigrip
	NobelReplace™ WP 5.0	35 Ncm	☆	Nobel Biocare® Unigrip
	NobelReplace™ 6.0	35 Ncm	☆	Nobel Biocare® Unigrip
	Brånemark System™ NP 3.3	35 Ncm	☆	Nobel Biocare® Unigrip
	Brånemark System™ RP 3.75, 4.0	35 Ncm	☆	Nobel Biocare® Unigrip
	Brånemark System™ WP 5.0, 5.5	35 Ncm	☆	Nobel Biocare® Unigrip
<b>Straumann®</b>				
	Bone Level 3.3 NC (Narrow CrossFit™)	35 Ncm	☆	SCS Screwdriver
	Bone Level 4.1, 4.8 RC (Regular CrossFit™)	35 Ncm	☆	SCS Screwdriver
	Standard/Standard Plus 4.8 RN (Regular Neck)	35 Ncm	☆	SCS Screwdriver
	Standard/Standard Plus 4.8 WN (Wide Neck)	35 Ncm	☆	SCS Screwdriver
	Standard Plus 3.3 NN (Narrow Neck)	35 Ncm	☆	SCS Screwdriver
<b>Zimmer Biomet</b>				
	3i Certain™ 3.25, 4.0, 4/3, 5.0, 5/4, 6.0	20 Ncm	⬡	0.048" Hex
	3i Certain™ Prevail™ 3/4/3, 4/5/4, 5/6/5, 4/3, 5/4, 6/5	20 Ncm	⬡	0.048" Hex
	3i Certain™ XP 4/5, 5/6	20 Ncm	⬡	0.048" Hex
	3i Miniplant™ 3.25	20 Ncm	⬡	0.048" Hex
	3i Osseotite™ 3.75, 4.0, 5.0, 6.0	20 Ncm	⬡	0.048" Hex
	3i Osseotite™ XP 4/5, 5/6	20 Ncm	⬡	0.048" Hex
	Zimmer Tapered Screw-Vent 3.5, 4.5, 5.7/Tapered Screw-Vent 3.5, 4.5, 5.7 friction fit	30 Ncm	⬡	1.25 mm Hex
	Zimmer Screw-Vent 3.3, 3.7, 4.5/Screw-Vent 3.3, 3.7, 4.5 friction fit	35 Ncm	⬡	1.25 mm Hex
	Zimmer SwissPlus™ 4.8	35 Ncm	⬡	1.25 mm Hex

\*These implant torque values are recommendations only and are provided for your reference. They are subject to change without notice in order to adhere to manufacturer's recommendations. No guarantee is implied or given for this information.

