

A Guide to Selecting Dentures



Simplicity Milled

Best for immediate placement or for budget conscious patients.



Integrity Line Milled

A balance of strength, fit, esthetics, and price.



Integrity Line Traditional

A classic solution balancing function, fit, and esthetics.



Ultra Line Traditional

Our most customized denture, offers the highest level of esthetics, fit, and function.



	<p>Simplicity Milled Best for immediate placement or for budget conscious patients.</p> 	<p>Integrity Line Milled A balance of strength, fit, esthetics, and price.</p> 	<p>Integrity Line Traditional A classic solution balancing function, fit, and esthetics.</p> 	<p>Ultra Line Traditional Our most customized denture, offers the highest level of esthetics, fit, and function.</p> 
Material and Finish	<p>Acrylic milled base MMA-free UV resin printed teeth Smooth finish Smooth palate</p>	<p>Ivotion® premium milled acrylic denture base Ivotion® Cross Linked PMMA milled teeth Naturally contoured gingiva Anatomical palate</p>	<p>Ivocap® acrylic Naturally contoured gingiva Anatomical palate</p>	<p>Ivocap® acrylic Naturally contoured gingiva Customized gingiva coloring available Anatomical palate</p>
Articulation	Digital	Digital	Stratos 100	Stratos 100
Set Up	Basic set up	Esthetic tooth arrangement	Esthetic tooth arrangement	Custom anterior esthetic tooth arrangement
Try In	Finished denture at try in	Monolithic A2 denture try in	Wax try in	Wax try in
Resets Included	No	1 Reset - if necessary	1 Reset - if necessary	2 Resets - if necessary
Base and Bite Registration	Printed base, bite block	Printed base, bite block	Processed or light cure base	Processed or light cure base
Personal ID Included	No	No	No	Yes
Tooth Shades Available	A1, A2, A3, A3.5, B1, Bleach	A1, A2, A3, A3.5, B1, Bleach	See Reverse for Tooth Options	See Reverse for Tooth Options
Warranty	1 Year	1 Year	1 Year	1 Year
Price	\$259 per arch (includes teeth)	\$425 per arch (includes teeth)	\$350 per arch (not including teeth)	\$450 per arch (not including teeth)



Traditional Tooth Options (Integrity Line Traditional and Ultra Line Traditional only)

	Ivostar®/Gnathostar® <i>For value-driven restorations.</i>	Vivodent® DCL <i>For any removable restorations.</i>	Phonares® II <i>For implant restorations or highly esthetic and high wear considerations.</i>
Material	Cross Linked PMMA	Double Cross Linked (DCL)	Nano Hybrid Composite (NHC)
Shade	20 Chromoscop shades 16 A-D shades, Bleach BL1 to BL4	16 A-D shades, Bleach BL1 to BL4	16 A-D shades, Bleach BL1 to BL4
Anterior Moulds	15 maxillary, 6 mandibular	42 maxillary, 8 mandibular	18 maxillary, 6 mandibular
Price (1x14)	\$64	\$148	\$186

ADA Codes:

Complete Denture - Maxillary D5110
Complete Denture - Mandibular D5120

Immediate Denture - Maxillary D5130
Immediate Denture - Mandibular D5140

Adjust Complete Denture - Maxillary D5410
Adjust Complete Denture - Mandibular D5411

Repair Broken Complete Denture Base - Maxillary D5512
Repair Broken Complete Denture Base - Mandibular D5511

Replace Missing or Broken Teeth - Complete Denture (each tooth) - D5520

Rebase Complete Denture - Maxillary D5710
Rebase Complete Denture - Mandibular D5711

Reline Complete Denture (direct) - Maxillary D5730
Reline Complete Denture (direct) - Mandibular D5731

Reline Complete Denture (indirect) - Maxillary D5750
Reline Complete Denture (indirect) - Mandibular D5751

Interim Complete Denture - Maxillary D5810
Interim Complete Denture - Mandibular D5811



Pricing effective 04/25/2025. Pricing subject to change.



Dental Prosthetic Services

1900 51st Street NE
Cedar Rapids, IA 52402
800-332-3341 319-393-1990
www.DPSdental.com

