

Time Schedule



Number of **business days** required for production.
(Does not include pick up or delivery time)

Crown & Bridge - Kami Waltemeyer

4 days	Nautilus Full-Contour (3 units or less)
4 days	Dura-Temp bridge
5 days	Nautilus Elite Multi
5 days	Gradia® Composite
6 days	Full cast metal
6 days	Nautilus Core
6 days	Nautilus Full-Contour (4 units or more)
6 days	Simple Z
6 days	All-ceramic, IPS e.max®
9 days	Porcelain fused to metal crown or bridge
Gallery Collection - Call Troy Carn for time schedule and case planning	

Removable - Jeanie Wood

1-2 days	Acrylic repairs/addition without metal
2-3 days	Repair/addition with metal
3 days	Soft reline

Traditional Dentures

2 days	Bite registrations
2 days	Impression trays
3 days	Set-up
4 days	Finish
4 days	Immediate dentures
Add one day - silicone soft liner/gasket	

Value & Integrity Milled Dentures

5 days	Bite registrations
2 days	Impression trays
5 days	Set-up
5 days	Finish

Traditional Partial

2 days	Design and estimate
3 days	Repairs
6 days	Cast frame only
8 days	Cast frame for tooth try-in
10 days	Cast frame and acrylic finish
5 days	Flexible partials
4 days	Flippers/acrylic partials
Add one day - silicone soft liner/gasket	

Simplicity Acrylic Partial

4 days	Simplicity Acrylic Partial
--------	----------------------------

Orthodontics & Sleep Appliances - Colleen Kesterson

2 days	Space maintaining appliances	Simply Printed Splints
2 days	Proform mouthguard/nightguard	4 days Simply Printed Splint
3 days	Hawley retainers	
5 days	Functional appliances	
3 days	Ortho repairs	
3 days	Bite splints	
5-7 days	Sleep appliance repairs	
10 days	Obstructive sleep apnea appliances	
15 days	OrthoApnea NOA and Panthera sleep apnea appliances	
14 days	Obstructive sleep apnea appliance demo models	

Implants - Mark Whittaker and Dennis Becker, CDT

Call for time schedule and case planning

Updated 04/21/2026

Product Criteria



Nautilus Core

Porcelain-fused-to-zirconia

- 1.5 mm occlusal reduction
- Any margin design
- Avoid sharp or right angles and undercuts
- Anterior crowns and bridges
- Conventional cementation

Nautilus Elite Multi

Solid multilayer zirconia

- 1.0 mm occlusal reduction
- Any margin design
- Avoid sharp or right angles and undercuts
- Anterior crowns and bridges
- Conventional cementation

Nautilus Full-Contour

Solid zirconia

- 1.0 mm occlusal reduction
- Any margin design
- Avoid sharp or right angles and undercuts
- Posterior crowns and bridges
- Conventional cementation

Simple Z

Solid zirconia

- 1.0 mm occlusal reduction
- Any margin design
- Avoid sharp or right angles and undercuts
- Posterior single-unit crowns only
- Available in Classic shades
- Conventional cementation

IPS e.max®

- 2.0 mm occlusal/axial reduction
- 1.0 mm (90°) shoulder margins
- Avoid sharp or right angles and undercuts
- Anterior crowns
- Conventional cementation or bonding

IPS e.max® Veneers

- 1.0 mm - 1.5 mm incisal reduction
- Correct preparation of the chamfer margins interproximally
- Depth cuts of 0.6mm - 0.8 mm
- Must be bonded

Gradia® Composite

Polyglass for inlay, onlay, and crowns

- 1.5 mm occlusal/axial reduction
- Avoid sharp or right angles and undercuts
- Must be bonded

DPS PFM

- 1.5 mm occlusal reduction
- Any margin design
- Any collar design
- Anterior and posterior crowns and bridges
- Conventional cementation



Dental Prosthetic Services, Inc
1900 51st Street NE
Cedar Rapids, Iowa 52402
319-393-1990 or 800-332-3341
Fax: 319-393-8455
www.DPSdental.com